


AI SIMÔ



AI SIMO CORPORATION CO., LTD



AISIMO autosampler vials

AISIMO autosampler vials improve your laboratory productivity by reducing costs and saving time. The cleanest on the market, these products considerably reduce the risk of test results compromised by ghost peaks, damaged needles or dislodged septa, decreasing analysis failures and sample reruns. AISIMO autosampler vials and caps are supplied certified as a standard at no extra cost, ensuring that with vials from AISIMO you will achieve consistent, reliable and reproducible results time after time.

All of our vials, inserts, caps, septa are extensively quality controlled and many are available as certified product. They are fully guaranteed to be reproducible and compatible with your autosampler.

We offer a full range of 2ml vials, 4ml vials and 10ml vials, up to 20ml vials and above, including VOA and EPA vials at 40mls. These are available in Glass Vial and also Plastic Vial format for many types.

Glass— clear glass vials are made from Type I, 33 expansion glass, with very low free ion content
Dimensional—vials are held to tight dimensional tolerances to avoid damage to the needle or instrument

Cap/Septa— tested by head space GC to ensure proper curing of silicone polymer

Vials, caps and septa – a wide range suitable for any application, and pre-assembled kits for your convenience.

AISIMO appreciate that selecting the correct chromatography vial is not always straight forward which is why our renowned customer service is here to help.

Order information:


Part No.	Products Name/Size	Unit/pk
2mm Snap neck vial		
11AVC-SN	2ml Clear vial, snap top	100
11AVC-SW	2ml Clear vial, snap top, graduated with write-on spot	100
11AVA-SN	2ml Amber vial, snap top	100
11AVA-SW	2ml Amber vial, snap top, graduated with write-on spot	100
Caps + Septa preassembled		
CS11-SW	White PTFE/red silicone septa + PE transparent cap with hole, for 2ml snap top vial	100
CS11-SR	Red PTFE/white silicone septa + PE transparent cap with hole, for 2ml snap top vial	100
8-425 screw neck vial		
8AVC-N	2ml Clear vial, 8-425 screw top	100

8AVC-W	2ml Clear vial, 8-425 screw top, graduated with write-on spot	100
8AVA-N	2ml Amber vial, 8-425 screw top	100
8AVA-W	2ml Amber vial, 8-425 screw top, graduated with write-on spot	100
Caps + Septa preassembled		
CS08-W	White PTFE/red silicone septa + Black screw cap with hole, for 2ml 8-425 screw top vial	100
CS08-R	Red PTFE/white silicone septa + Black screw cap with hole, for 2ml 8-425 screw top vial	100
9mm screw neck vial		
9AVC-N	2ml Clear vial, 9mm screw top	100
9AVC-W	2ml Clear vial, 9mm screw top, graduated with write-on spot	100
9AVA-N	2ml Amber vial, 9mm screw top	100
9AVA-W	2ml Amber vial, 9mm screw top, graduated with write-on spot	100
Caps + Septa preassembled		
CS09-W	White PTFE/red silicone septa + Blue screw cap with hole, for 2ml 9mm screw top vial	100
CS09-R	Red PTFE/white silicone septa + Blue screw cap with hole, for 2ml 9mm screw top vial	100
CS09-P	Blue PTFE/white silicone septa, Pre-slit + Blue screw cap with hole, for 2ml 9mm screw top vial	100
10-425 screw neck vial		
10AVC-W	2ml Clear vial, 10-425 screw top, graduated with write-on spot	100
10AVA-W	2ml Amber vial, 10-425 screw top, graduated with write-on spot	100
Caps + Septa preassembled		
CS10-W	White PTFE/red silicone septa + Black screw cap with hole, for 2ml 10-425 screw top vial	100
CS10-R	Red PTFE/white silicone septa + Black screw cap with hole, for 2ml 10-425 screw top vial	100
11mm Crimp neck vial		
11AVC-CN	2ml Clear vial, crimp top	100
11AVC-CW	2ml Clear vial, crimp top, graduated with write-on spot	100
11AVA-CN	2ml Amber vial, crimp top	100
11AVA-CW	2ml Amber vial, crimp top, graduated with write-on spot	100
Caps+ Septa preassembled		

CS11-CW	White PTFE/red silicone septa+ Aluminium cap, for 2ml crimp top vial	100
CA11-CR	Red PTFE/white silicone septa+ Aluminium cap, for 2ml crimp top vial	100
4ml 13-425 Screw neck vial		
13AVC-CN	4ml Clear vial, 13-425 screw top	100
13AVC-CW	4ml Clear vial, 13-425 screw top, graduated with write-on spot	100
13AVA-CN	4ml Amber vial, 13-425 screw top	100
13AVA-CW	4ml Amber vial, 13-425 screw top, graduated with write-on spot	100
Caps + Septa preassembled		
CS13-CW	White PTFE/red silicone septa + Black screw cap with hole, for 4ml 13-425 screw top vial	100
CS13-CR	Red PTFE/white silicone septa + Black screw cap with hole, for 4ml 13-425 screw top vial	100
18mm Screw neck Headspace vial		
18GC46CF	10ml Clear vial, 22.5×46mm, 18mm screw top, flat bottom	100
18GC75CF	20ml Clear vial, 22.5×75mm, 18mm screw top, flat bottom	100
18GC46AF	10ml Amber vial, 22.5×46mm, 18mm screw top, flat bottom	100
18GC75AF	20ml Amber vial, 22.5×75mm, 18mm screw top, flat bottom	100
Caps + Septa preassembled		
CS18MM	Molded PTFE/butyl septa + Magnetic cap with hole, for 18mm screw top vial	100
CS18WM	White PTFE/white silicone septa + Magnetic cap with hole, for 18mm screw top vial	100
20mm Crimp neck Headspace vial		
20GC38CF	6 ml Clear vial, 22.5 x 38mm, crimp top, flat bottom	100
20GC46CF	10ml Clear vial, 22.5×46mm, crimp top, flat bottom	100
20GC 75CF	20ml Clear vial, 22.5×75mm, crimp top, flat bottom	100
20GC 38CR	6 ml Clear vial, 22.5 x 38mm, crimp top, rounded bottom	100
20GC 46CR	10ml Clear vial, 22.5×46mm, crimp top, rounded bottom	100
20GC 75CR	20ml Clear vial, 22.5×75mm, crimp top, rounded bottom	100

20GC46AR	10ml Amber vial, 22.5×46mm, crimp top, rounded bottom	100
20GC75AR	20ml Amber vial, 22.5×75mm, crimp top, rounded bottom	100
Caps + Septa preassembled		
CS20SL	Stopper, Molded Gray Butyl, Locking Design, for 20mm crimp vial	100
CS20WA	White PTFE/white silicone septa + Aluminium cap, for 20mm crimp vial	100
CS20RA	Red PTFE/white silicone septa + Aluminium cap, for 20mm crimp vial	100
CS20MA	Molded PTFE/butyl septa + Aluminium cap, for 20mm crimp vial	100

Other specifications are available on request.



More service and products:

info@aisimo.com

www.aisimo.com

+44-20-32871993 (England)

+86-21-31131618 (China)

+1-518-8939698 (America)

+61-28-0036806 (Australia)

AIMÔ

To place an order, please send item number and quantity to info@aisimo.com